





BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-1270-50F	mm	70	12.5	1.27	1.27
	in	2.76	0.49	0.050	0.050
Replaces N/A					
					<p>No Tooth Offset</p> <p>Flared Blade</p>
K5-1290-47	mm	90	12.5	1.19	1.19
	in	3.54	0.49	0.047	0.047
Replaces 05.003.1045					
					No Tooth Offset
K5-1290-50	mm	90	12.5	1.27	1.27
	in	3.54	0.49	0.050	0.050
Replaces 05.003.1065					
					No Tooth Offset
K5-1290-54	mm	90	12.5	1.37	1.37
	in	3.54	0.49	0.054	0.054
Replaces 05.003.1075					
					No Tooth Offset

Compatible with DePuy® Synthes®
Large Bone Oscillating Saw System



BUSA®
Proprietary

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Compatible with DePuy® Synthes®
Large Bone Oscillating Saw System







BUSA®
Proprietary

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-1990-39	mm	90	19	1.00	1.00
	in	3.54	0.75	0.039	0.039
Replaces 05.003.1125					
					No Tooth Offset
K5-1990-47	mm	90	19	1.19	1.19
	in	3.54	0.75	0.047	0.047
Replaces 05.003.1145					
					No Tooth Offset
K5-1990-50	mm	90	19	1.27	1.27
	in	3.54	0.75	0.050	0.050
Replaces 05.003.1165					
					No Tooth Offset
K5-1990-54	mm	90	19	1.37	1.37
	in	3.54	0.75	0.054	0.054
Replaces 05.003.1175					
					No Tooth Offset

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-2190-39F	mm	90	21	0.99	0.99
	in	3.54	0.83	0.039	0.039
Replaces 05.002.3065					
					No Tooth Offset Flared Blade
K5-2190-47F	mm	90	21	1.19	1.19
	in	3.54	0.83	0.047	0.047
Replaces 05.002.3025					
					No Tooth Offset Flared Blade
K5-2190-50F	mm	90	21	1.27	1.27
	in	3.54	0.83	0.050	0.050
Replaces 05.002.3035					
					No Tooth Offset Flared Blade
K5-2190-54F	mm	90	21	1.37	1.37
	in	3.54	0.83	0.054	0.054
Replaces 05.002.3045					
					No Tooth Offset Flared Blade

Compatible with DePuy® Synthes®
Large Bone Oscillating Saw System



BUSA®
Proprietary

BUSA® Large Bone Saw Blades - Oscillating BR2000™ (Aggressive Cut) - Proprietary

Compatible with DePuy® Synthes® Large Bone Oscillating Saw System

Compatible with DePuy® Synthes®
Large Bone Oscillating Saw System



BUSA®
Proprietary

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-2590-39	mm	90	25	1.00	1.00
	in	3.54	0.98	0.039	0.039
Replaces 05.003.1325					
					No Tooth Offset Flared Blade
K5-2590-47	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
Replaces 05.003.1345					
					No Tooth Offset
K5-2590-50	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.050
Replaces 05.003.1365					
					No Tooth Offset
K5-2590-54	mm	90	25	1.37	1.37
	in	3.54	0.98	0.054	0.054
Replaces 05.003.1375					
					No Tooth Offset