





BUSA® Large Bone Saw Blades - Sagittal BR2000™ (Aggressive Cut) - Proprietary

Compatible with Stryker® Sagittal Saw Systems: 7, 6, 5, 4, 2000 & EHD





Compatible with Stryker® Sagittal Saw Systems:
7, 6, 5, 4, 2000 & EHD

BUSA®
Proprietary

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-1290-35	mm	90	12.5	0.89	0.89
	in	3.54	0.49	0.035	0.035
<i>Substitute for Stryker® 6113-089-090</i>					
					No Tooth Offset
BR1-1290-39	mm	90	12.5	0.99	0.99
	in	3.54	0.49	0.039	0.039
<i>Substitute for Stryker® 6113-097-090</i>					
					No Tooth Offset
BR1-1290-47	mm	90	12.5	1.19	1.19
	in	3.54	0.49	0.047	0.047
<i>Substitute for Stryker® 6113-119-090 & 2108-352</i>					
					No Tooth Offset
BR1-1290-50	mm	90	12.5	1.27	1.27
	in	3.54	0.49	0.050	0.050
<i>Substitute for Stryker® 6113-127-090</i>					
					No Tooth Offset

BUSA® Large Bone Saw Blades - Sagittal BR2000™ (Aggressive Cut) - Proprietary

Compatible with Stryker® Sagittal Saw Systems: 7, 6, 5, 4, 2000 & EHD

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-1290-54	mm	90	12.5	1.37	1.37
	in	3.54	0.49	0.054	0.054
<i>Substitute for Stryker® 6113-137-090</i>					
					No Tooth Offset
BR1-1990-35	mm	90	19	0.89	0.89
	in	3.54	0.75	0.035	0.035
<i>Substitute for Stryker® 6118-089-090</i>					
					No Tooth Offset
BR1-1990-39	mm	90	19	0.99	0.99
	in	3.54	0.75	0.039	0.039
<i>Substitute for Stryker® 6118-097-090</i>					
					No Tooth Offset
BR1-1990-47	mm	90	19	1.19	1.19
	in	3.54	0.75	0.047	0.047
<i>Substitute for Stryker® 6118-119-090</i>					
					No Tooth Offset

Compatible with Stryker® Sagittal Saw Systems:
7, 6, 5, 4, 2000 & EHD







BUSA®
Proprietary

BUSA® Large Bone Saw Blades - Sagittal BR2000™ (Aggressive Cut) - Proprietary

Compatible with Stryker® Sagittal Saw Systems: 7, 6, 5, 4, 2000 & EHD





Compatible with Stryker® Sagittal Saw Systems:
7, 6, 5, 4, 2000 & EHD

BUSA®
Proprietary

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-19100-47	mm	100	19	1.19	1.19
	in	3.94	0.75	0.047	0.047
<i>Substitute for Stryker® 6118-119-100</i>					
 <p style="text-align: right;">No Tooth Offset</p>					
BR1-1990-50	mm	90	19	1.27	1.27
	in	3.54	0.75	0.050	0.050
<i>Substitute for Stryker® 6118-127-090 & 2108-382</i>					
 <p style="text-align: right;">No Tooth Offset</p>					
BR1-19100-50	mm	100	19	1.27	1.27
	in	3.94	0.75	0.050	0.050
<i>Substitute for Stryker® 6118-127-100 & 2108-385</i>					
 <p style="text-align: right;">No Tooth Offset</p>					
BR1-19110-50	mm	110	19	1.27	1.27
	in	4.33	0.75	0.050	0.050
<i>Substitute for Stryker® 6118-127-110</i>					
 <p style="text-align: right;">No Tooth Offset</p>					

BUSA® Large Bone Saw Blades - Sagittal BR2000™ (Aggressive Cut) - Proprietary

Compatible with Stryker® Sagittal Saw Systems: 7, 6, 5, 4, 2000 & EHD

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-1990-54	mm	90	19	1.37	1.37
	in	3.54	0.75	0.054	0.054
<i>Substitute for Stryker® 6118-137-090</i>					
					No Tooth Offset
BR1-19100-54	mm	100	19	1.37	1.37
	in	3.94	0.75	0.054	0.054
<i>Substitute for Stryker® 6118-137-100</i>					
					No Tooth Offset
BR1-2190-39F	mm	90	21	0.99	0.99
	in	3.54	0.83	0.039	0.039
<i>Substitute for Stryker® 6221-100-090</i>					
					No Tooth Offset Flared Blade
BR1-2190-47F	mm	90	21	1.19	1.19
	in	3.54	0.83	0.047	0.047
<i>Substitute for Stryker® to 6221-119-090</i>					
					No Tooth Offset Flared Blade

Compatible with Stryker® Sagittal Saw Systems:
7, 6, 5, 4, 2000 & EHD







BUSA®
Proprietary

BUSA® Large Bone Saw Blades - Sagittal BR2000™ (Aggressive Cut) - Proprietary

Compatible with Stryker® Sagittal Saw Systems:
7, 6, 5, 4, 2000 & EHD





BUSA®
Proprietary

Compatible with Stryker® Sagittal Saw Systems: 7, 6, 5, 4, 2000 & EHD

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-2190-50F	mm	90	21	1.27	1.27
	in	3.54	0.83	0.050	0.050
<i>Substitute for Stryker® 6221-127-090, 2108-382-001 & 2108-382-3</i>					
					No Tooth Offset Flared Blade
BR1-2190-54F	mm	90	21	1.37	1.37
	in	3.54	0.83	0.054	0.054
<i>Substitute for Stryker® 6221-137-090 & 2108-382-2</i>					
					No Tooth Offset Flared Blade
BR1-2590-35	mm	90	25	0.89	0.89
	in	3.54	0.98	0.035	0.035
<i>Substitute for Stryker® 6125-089-090 & 2108-302</i>					
					No Tooth Offset
BR1-2590-39	mm	90	25	0.99	0.99
	in	3.54	0.98	0.039	0.039
<i>Substitute for Stryker® 6125-097-090 & 2108-393-006</i>					
					No Tooth Offset

BUSA® Large Bone Saw Blades - Sagittal BR2000™ (Aggressive Cut) - Proprietary

Compatible with Stryker® Sagittal Saw Systems: 7, 6, 5, 4, 2000 & EHD

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-2590-39F	mm	90	25	0.99	0.99
	in	3.54	0.98	0.039	0.039
<i>Substitute for Stryker® 2108-393-005</i>					
					<p>No Tooth Offset</p> <p>Flared Blade</p>
BR1-2590-47	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
<i>Substitute for Stryker® 6125-119-090 & 2108-351</i>					
					No Tooth Offset
BR1-2590-47F	mm	90	25	1.19	1.19
	in	3.54	0.98	0.047	0.047
<i>Substitute for Stryker® 2108-351-002</i>					
					<p>No Tooth Offset</p> <p>Flared Blade</p>
BR1-2590-50	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.050
<i>Substitute for Stryker® 6125-127-090, 2108-376-002 & 2108-377</i>					
					No Tooth Offset

Compatible with Stryker® Sagittal Saw Systems:
7, 6, 5, 4, 2000 & EHD



BUSA®
Proprietary

BUSA® Large Bone Saw Blades - Sagittal BR2000™ (Aggressive Cut) - Proprietary

Compatible with Stryker® Sagittal Saw Systems: 7, 6, 5, 4, 2000 & EHD

Compatible with Stryker® Sagittal Saw Systems:
7, 6, 5, 4, 2000 & EHD






BUSA®
Proprietary

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-25100-50	mm	100	25	1.27	1.27
	in	3.94	0.98	0.050	0.050
Substitute for Stryker® 6125-127-100 & 2108-328					
BR1-2590-50F	mm	90	25	1.27	1.27
	in	3.54	0.98	0.050	0.050
Substitute for Stryker® 2108-383					
BR1-2590-54	mm	90	25	1.37	1.37
	in	3.54	0.98	0.054	0.054
Substitute for Stryker® 6125-137-090					
BR1-2590-54F	mm	90	25	1.37	1.37
	in	3.54	0.98	0.054	0.054
Substitute for N/A					

BUSA® Large Bone Saw Blades - Sagittal BR2000™ (Aggressive Cut) - Proprietary

Compatible with Stryker® Sagittal Saw Systems: 7, 6, 5, 4, 2000 & EHD

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
BR1-2590-58	mm	90	25	1.47	1.47
	in	3.54	0.98	0.058	0.058
<i>Substitute for Stryker® 6125-147-090</i>					
					No Tooth Offset
BR1-2590-58F	mm	90	25	1.47	1.47
	in	3.54	0.98	0.058	0.058
<i>Substitute for N/A</i>					
					No Tooth Offset Flared Blade
BR1-3165-39F	mm	65	31	0.99	0.99
	in	2.56	1.22	0.039	0.039
<i>Substitute for N/A</i>					
					No Tooth Offset Flared Blade

Compatible with Stryker® Sagittal Saw Systems:
7, 6, 5, 4, 2000 & EHD



BUSA®
Proprietary